



Leading homeopathic software
radaropus
THE WORLD'S No.1 SINCE 1982

RadarOpus 2.1

Installation

on Windows or Mac

**With update information
from RadarClassic to RadarOpus**

Table of content

RADAROPUS NEW INSTALLATION	3
New Installation via Download	3
Next installation steps	3
USER ACCOUNT	6
User Account (Local & RadarOpus Cloud)	6
Learn to use RadarOpus	7
Advised program settings	7
How to move your SOFT-KEY license to another computer	7
Technical support	8
Import patient information from an old RadarClassic	8
Installation via a Memory Stick	8
Installation of the Live update 2.1 on top of Opus 2.0.xx	9
How to transfer a SOFT-KEY license to a different computer	9
How to transfer a HARD-KEY license to a different computer	9
TECHNICAL SUPPORT	11
Technical support & solutions	11
Compatibility with High-Resolution screens	13
How to contact technical support	15
RE-INSTALLATION RADAROPUS	16
Re-Install RadarOpus	16
UPDATE RADARCLASSIC TO RADAROPUS	17
Update from RadarClassic to RadarOpus	17
Import Patients, Clipboards, Additions from RadarClassic	17

RadarOpus New Installation

New Installation via Download



- ⇒ If your internet speed is not fast enough to download the installation file, please continue to read the chapter: **Installation via a Memory Stick**.
- ⇒ If you want to Update or do a Re-installation, please read the chapter: **Re-installation and Update instructions**.

New installation:

- Go to www.RadarOpus.com
- Click in the menu bar on: **DOWNLOAD**
- Then at the bottom, select:
DOWNLOAD RADAROPUS
- Then select: **DOWNLOAD FULL**
- Select the correct link for:
 - **Mac**
 - **Windows 8.1 / 10**
 - **Windows 7**
- Save the file in a location you can find back easily, for example in the **folder Downloads**
Wait until the download is ready.
- Then, if you are running an **anti-virus** then it is best to **turn this off temporarily**.
- Then open the folder **downloads**:


On Mac: **Double click** on the downloaded file: **RadarOpus-OSX-installer.dmg**
Then you will see several folders and files, then **double click** on the installer file

On Windows: **Right mouse click** on the downloaded file: **RadarOpus_2.1.xx_windows-installer.iso**
Then select **Open with**
Then select **Explorer**
You will see several folders and files, then **double click** on the installer file




Next installation steps

- Select the **Language**, then click **Next**
(This is the language of the menu bar, it is not related to the language of the repertory)




- Accept the **License agreement** and then click on the button **Next**



- A. Select the **Drive letter** where you want to install RadarOpus (only possible on Windows).

(Select a Drive letter where you have minimal 30 GB free)

- B. Click **Next** to start the installation. The installation can take from 10 minutes up to 45 minutes, depending on the speed of the computer and if you have an SSD or an older HDD.



- C. After the installation is ready, **restart the computer**.

- D. Start your RadarOpus via the program icon



- E. You will see a window with two choices. Click on the correct box: did your package come with a:

HARD-KEY (usb-key) protection

Or:

SOFT-KEY (CmActLicense) protection




If you have a HARD-KEY, plug this **mini USB key** in the computer (also called “Wibu Code meter”)

- It is best to plug in the usb key directly into your computer, not via a usb-hub, or cable.
- Wait a few seconds, before you continue.

F. Registration

Enter the following information

- **First + Last name**
 - **Email address**,
 - Select your **country**
 - Your RadarOpus **license number**
 - Your **Activation key**
- (Start at the beginning of the field to avoid a space)



- G. Click **Next**

If a message appears **“Cannot connect to the server”** then **please wait 60 seconds**.

The install program will try several times, and it may take 60 seconds before the connection is established. A “Thumbs Up” icon will appear when ready.

The installation is now complete, and you can **start** RadarOpus with the desktop icon:



If it still does not work, please contact us for support. Do not install again.

Please read the section: ***How to contact technical support***

H. **Starting RadarOpus: Create a user account**

In many countries there are strict rules regarding patient information on a computer (GDPR, AVG, or HIPAA).

RadarOpus complies with these rules if you create an account.

Please read below for more information

Attention: If you want to move your software license to a different computer, please read:

How to move your SOFT-KEY license to another computer

(soon a new website with the license manager will be available)

I. Updates

Probably you will receive immediately a message that there is a update available e.g. 2.1.12

You can install it.

There are also Content updates available.

You can install all the Free items.

Attention: the Repertory ADDITIONS files, each will take about 10 min to install.

User Account

User Account (Local & RadarOpus Cloud)


Step 1:

After installing the update, and starting RadarOpus you will be asked if you want to: **create an account?**

It is **not yet mandatory** to create an account.

If you are a professional practitioner we advise you to do so because in many countries new regulations oblige you to protect all software which contains patient information with a **password** and **database encryption** (HIPAA regulations, GDPR, AVG).

You are not allowed anymore to have e.g. Word files with confidential patient information in them stored on your computer, nor is it allowed to store such files in a cloud drive.



30 days reminder:

You can select: **No, remind me later.**

For 30 days you will receive a reminder to create an account.

If you have not done so after 30 days, you will be able to continue to work without a password.



We are studying the regulations, and maybe in future it will become mandatory to create an account!

It is possible to create an account later at any moment: In the RadarOpus manual you can read more how to do this.

Step 2: Creating an account

If you select: **Yes, create an account now**, there are two parts:

Local account

→ You can create only a **Local account** to protect your RadarOpus patient information.

A "local account" is just the **login** when you start RadarOpus.


- Please enter: **First & Last name**, and a **Password**
- Select a security question & give an answer

RadarOpus cloud account (Clificol)

→ You can also select to create a RadarOpus-Cloud account.

- You will receive **Clificol access** (this is a cured cases database). For more information see www.clificol.net

- In future versions of RadarOpus more functions will come available.



Attention: Check the Account.txt file


On your **desktop** will be placed a file called **Account.txt**

It will contain the details of the two accounts you created:

- The password for the **Local account**
- A password for the **Clificol (RadarOpus Cloud)**
(this is different from your local account password)

Print or note down the login details.

And make sure you **delete** this Account.txt file, and also delete it from your trash box.



Confirm the Clificol account

You will receive an **email** asking you to confirm you RadarOpus Cloud (Clificol) account.

Use the password mentioned in the **Account.txt** document

Learn to use RadarOpus

- A. You can start by viewing the: **Introduction Tutorial Videos**.
Ask you representative for the correct link, you can find many in the Academy on www.radaropus.com
- B. There is an additional **RadarOpus Manual** or consult the **Help on line**
- C. Make sure you check out the document **Advised program settings** for the most optimal settings for your country/language.

Advised program settings

Make sure you ask your local representative for the **advised program settings** which will work best for your country and language.

How to move your SOFT-KEY license to another computer

From RadarOpus version 2.2 you can easily move your license to another (new) computer via the LicenseServer-website.

With version 2.1. you need to ask you representative for instructions.

Attention

*After you installed RadarOpus with **SOFT-KEY** you can only move it to another computer after you de-activate it from the first computer via a website (license server).*

De-activate it for example before you sell or throw away your computer with a RadarOpus Soft-key license installed on it, or before you give this computer for repair.

If you have a **HARD-key** you can simply unplug the usb key and place it on the new computer and install the full Installer on that new computer.



Technical support

If there is any problem with the installation of registration, please continue to read the chapter *Technical support & solutions*

Import patient information from an old RadarClassic

If you need to import **Cases** or **Winchip Patient information** from RadarClassic 10.x to RadarOpus, please read the chapter: *Import Patients, Clipboards, Additions from RadarClassic*.

If you need help please contact your local RadarOpus representative or the Technical support.

Installation via a Memory Stick

If your internet speed is not fast enough to download the RadarOpus installation file, then please contact your local representative to order a memory stick.

- Place the memory stick in a usb port.
- Then continue to read the normal installation instructions

Attention if there is a problem

If you encounter a problem during installation from a memory stick, then first copy all files to the harddrive (for example to the folder \Downloads), then follow the above INSTALL instructions.

Summary of installation Windows or Mac

- Step 1: **Download** and install RadarOpus from: www.radaropus.com
(if your internet is too slow you can receive a memory stick to install from)
- Then continue to read the normal installation instructions

Re-installation and Update instructions

Please continue to read the chapter about: *Re-installation and Update instructions*

It is always best to download the most recent version of RadarOpus. Do not use old program versions anymore which you might still have on a DVD's or memory stick.

If you might still have a very old **Radar** 10 or 10.5 you do not need to de-install this first. The RadarOpus will be installed separately.

Installation of the Live update 2.1 on top of RadarOpus 2.0.xx

If you installed a RadarOpus version 2.0.xx your RadarOpus will give you automatically a message to install a live update version 2.1.11

It is also possible to download the update file manually from www.radaropus.com.

After download double click to start the installation.

How to transfer a SOFT-KEY license to a different computer

- Make a RadarOpus **backup** of your patient information, and save it to a location you can find back easily e.g. Documents or to an external memory stick or external Harddisk.
- (coming with RadarOpus version 2.2). Go to the **license server** (website link coming) and **de-activate** your license on this current computer.
- **Install RadarOpus** on the new computer according to the above installation instructions.
- Start RadarOpus, but do not create an account.
- **Import your RadarOpus backup** into this newly installed program:

*From the menu bar **Options**, then **Database management**, **Restore** .*

- Restart RadarOpus, and you need to use your previous used **Login** and **Password**.
(the login and password is stored inside the database)


How to transfer a HARD-KEY license to a different computer

If you have a HARD-KEY RadarOpus license with a user account created.

You can transfer your “User Account” and patient information to a new computer as follows:

- **If you have already saved patient info (and/or additions)**
If you had already saved patient information and/or personal additions on the previous computer, then make a **RadarOpus backup** on that computer, and **restore** it in the new computer, as follows:

- Make a RadarOpus **Backup** on the old computer.
- **Install RadarOpus** on the new computer according to the above installation instructions.
- Start RadarOpus but do not create an account. Select: **No, remind me later**



- Via the menu bar **Options**, **Database management**, **Restore** your RadarOpus backup. Then restart RadarOpus.
You then must login with your previous Login name and Password.


For more information, how to make a backup, and how to restore it, please see below.

- **If you had not yet saved any patient information** (or any personal additions).

Install RadarOpus on the new computer according to the above installation instructions.

You can create a (new) account and use the same (or any new) Login name and password if you like.

Select: **Yes, create an account**




You do **not** need to activate your Clificol (RadarOpus cloud) again.

Where can you find your old Login/Password details?

All the above **account information** you can find in the document **“Account.txt”** which you probably had printed and stored in a safe place, after you installed RadarOpus the first time on your previous computer.

You need to delete this file, and empty the trash, to be secure!



Technical support

Technical support & solutions

Error message: No License Found with a HARD-KEY

- If you start RadarOpus the first time, and you receive an error message: **No license found**

Step 1: Check if the RadarOpus USB* key (Code Meter) is placed (directly) on the computer, then do the following: (* This is not the same key as you used for the old Radar)

Step 2:

On Windows:


- Turn off the anti-virus program and/or the firewall to allow registration to contact internet.
- Open the Windows Explorer, Select **This PC**, then **Local Disk (C:)**, then open the folder **RadarOpus** and double click on the file called Activator.

On Mac:

- Via the **Apple menu**, go to System preferences. Then select **Security & Privacy**.

Change the security settings to: **“Allow apps downloaded from: Anywhere”**.

If you still get a warning message, select **“Open anyway”**



- Then go with the Finder to the **Program folder** and open the folder RadarOpus. Then double click on the files called **Activator**.

Step 3:

- If the above has not solved the problem, then **update the driver software** for the Wibu CodeMeter USB Key. *Please read the instructions below.*

Checklist for Windows users

- Make sure you have installed all the available **Windows updates**.
- Switch off temporarily the **Anti-Virus software** before you start the installation.
- If you have a **free 30 days Anti-Virus** software which period has expired, it often is necessary to fully uninstall this software and restart the computer. Sometimes this is also needed for your regular virus program, then after installation of RadarOpus you can install your regular virus scanner again or use the standard anti-virus from Windows itself which will switch on automatically.
- If you downloaded the Full RadarOpus installer as an ISO file, then do not double click on it, but use a **Right mouse click**, and select, **Open with**, then select **Windows explorer**


Check list for Mac users

- **MAC:** in rare situations, the Downloaded **ISO** file will not work. Then Download the **RadarOpus Zip file** and try to install this one. Then **Double click** on the Zip file after it has been downloaded. Then double click on the **RadarOpus OSX installer file**.
- It might happen that your Mac does not allow software to be installed which did not come from Apples App shop. Then follow the below steps:

Click on the **Apple menu**, then select **System Preferences**
Then select **Security & Privacy** icon.

Click on the **General** tab, then change in the section **“Allow apps downloaded from: Anywhere”**:
or if available click on: **“Open Anyway”**

Then **restart** your Mac and **plug in** the RadarOpus usb key and try to start the installation again.



Mac High Sierra, Mojave

- RadarOpus 2.1.11 is fully **High Sierra** and **Mojave** compatible.

Memory Stick issue

If you install from a memory stick, sometimes it is needed to first copy the content onto your local disk, before you double click on the installer file.

Installing the Wibu Code the Meter software and firmware update

Sometimes is it needed to manually download and install the most recent Wibu CodeMeter software.

Step 1: Download the Wibu software:

- Go to the website: **www.wibu.com**
- Then select: **Support and Downloads**, then **User**, then: **User software**
- Choose the download according to your operating system **Windows** or **Mac**.
- Open the Download folder, and **double click** on the downloaded file to install it.



Only for Hardware Key's:

Step 2: Update the USB key **firmware**:

- On **Mac**: go to the folder **Applications** and double click on **CodeMeter**
- On **Windows**: go to **START**, then **All programs**, then **Codemeter**, click on **Codemeter control center**


Or click on the CodeMeter icon in the system tray (bottom right corner of screen



Step 3:

- Click on the button with the arrow (open circle).

You will be told if there is a new version available or not. Click yes to update the firmware.




If you cannot find a solution, then please contact your local representative or the **Technical support**. See the chapter **How to contact the technical support**

Compatibility with High-Resolution screens

Radar Opus 2.1 has improved functionality to work with high resolution screens.


If you feel you like to improve the layout of RadarOpus on your computer, you can try the following to see which one will work optimal for your computer:

- Open in RadarOpus the **General settings**, and check ON or OFF “**Enable support for very high screen resolutions**” (tryout what works best).




- **On Mac** you can also go via the **Apple icon**, to **System Preferences, Displays**.

Change the **Default** for display to **Scaled** and tryout one higher or lower resolution setting.




On Windows Right mouse click on the desktop, then select **Display settings**.



Then first try to change “**Change the size of text**”, for example from 150% to **125%**, or **100%**.

Restart RadarOpus and see the result




If you cannot set the Font size and Windows sizes in a way you can comfortably work with RadarOpus, then please contact your local representative or the **Technical support**, or see the chapter **How to contact the technical support**

How to contact technical support

The fastest way to get help is to enter a technical support request via our **Ticket system**:

- Go to: www.radaropus.com
- Then select at the bottom of the page: **New support ticket**
- On working days, you will receive an email from technical support within 24 hours.
- The Technical support will often propose that you join them for a **support session**:
 - Go to the website support.zeus-soft.com
 - Click on the **name of the support person** who has contacted you.
 - Download the Bomgar software; if you cannot see the file download on your web page, go to your downloads folder (or wherever your computer downloads files from Internet) and double click on the Bomgar file (the icon is an orange square with a white B) to run it.
 - Click on Full Terms and Conditions and then on Accept
 - Click on Yes to accept the request for elevation of privileges
 - Click on Allow shared control of your computer



The technical support can only help you if you have a **Ticket** number

Re-Installation RadarOpus

Re-Install RadarOpus

If you have an **older** version of **RadarOpus** on your computer:

- **If it is an older version 1.40 or older:** It is best that everything of the previous version is **removed** before the new RadarOpus is installed.
 - Please contact the **support** who will do this for you.
 - Or you can do it yourself, read the document: [Clean installation](#)
- **If you have a version 1.41.x** on your computer running: you could only need to install the **live update 2.1** You can find this update on **www.radaropus.com**. But also best would be to do a full clean new installation.
- **If you already have a recent version 1.43.1 - 1.43.7** running: you may download and install the **live update 2.1** over your existing version 1.43.x. You can find this update on **www.radaropus.com**

New computer (or a full new installation):

- For any new installation do not use any older program DVD's. Always use the latest Full Installation file e.g. 2.1.x

Please read the: New installation instructions

Please also read the chapters above in this manual:

How to transfer a SOFT-KEY license to a different computer

How to transfer a HARD-KEY license to a different computer

Update RadarClassic to RadarOpus

Update from RadarClassic to RadarOpus

If you have important personal information stored in RadarClassic, it is advisable to make sure you **first make a copy of all personal information before you continue**.

And it will work easy if you keep the old version on your computer until RadarOpus is fully installed and your information

Please see information how to make a backup of RadarClassic.

You may leave your RadarClassic on your Windows computer when you install RadarOpus or remove it first. They are installed side by side and have no influence on each other.

RadarClassic normally can be found on the hard disk in the folder C:\Radar.

RadarOpus will be installed in the folder C:\RadarOpus and a second folder called C:\RadarOpus-PostgreSQL

If you have important personal information stored in RadarClassic, it is advisable to make sure you **first make a copy of all personal information before you continue**. Please see information how to make a backup of RadarClassic.

If you need help, please contact your local representative or contact the technical support

Import Patients, Clipboards, Additions from RadarClassic

Import (Transfer) information from Radar (8, 9, 10, 10.5) to RadarOpus

The following personal information in Radar and/or Winchip version can be converted to RadarOpus.

See video: [Import Cases](#) (RadarOpus v 1.30 of later)

- **Radar Clipboards (cases) – Catalogue*** (clipboards saved in the default way)
- **Radar Clipboards (cases) - Single Cas files*** (all clipboards saved as separate files on the HD)
- **Repertory additions****
- **Personal Keynotes, Freenotes, Personal EH-documents****
- **Winchip patient file***** (including optionally saved clipboards in Winchip)

(* **Note:** [Import on Apple Mac computers](#))

Your RadarOpus distributor will help with the conversion from the Windows Radar (8 -10.5) cases files (catalogue or single cases) to the Mac RadarOpus version. Only saved CAS files can be imported directly to RadarOpus-Mac. So, you will need to save your old Radar Catalogue of cases first to single **Cas** files.

If you have RadarOpus on Windows, you can import directly your Radar Catalogue of cases (technically 4 files). **Tip:** it is also possible to first import a Radar catalogue of cases into a **Windows-RadarOpus**, then make an Opus-backup - this Opus backup can be imported into a MAC-RadarOpus.

(**) **Note:** The import (transfer) of the Radar 9 or 10.x Personal Repertory additions is not yet possible in RadarOpus 1.43 or 2.0. Also, any personal Keynotes cannot yet be converted.

Please keep a backup copy of all your Radar information on a secondary storage media. Read the separate instruction: “How to make a backup of all your Radar personal information”.

You can find more information in the Radar 10.x Freenotes “The automatic backup” or see the “Repertory utilities manual”. RadarOpus version 2.1 will allow to convert your personal additions

Import Radar Clipboards

The clipboards (cases) in the old Radar program were by default stored in a so-called **Catalogue** (on the hard disk 4 files). It was also possible to save the cases as individual files, so-called **Cas-files**. For more information see below.

More explanation about “Catalogue” and “Cas”

The default way in the old Radar was to save the clipboards (cases) in a so-called ‘Catalogue’. The default catalogue name is: **Repert**. On the hard disk such a catalogue is located in the folder **c:\radar\cases** and consists of 4 files: **repert.da1, repert.da2, repert.id1, repert.id2**.

In Radar it was possible to save clipboards (cases) as individual files on the hard disk. They were usually stored in the folder **c:\radar\cases** and then in another subfolder, for example the folder **\Repert**. All the individual Cas files are located here **c:\radar\cases\repert** (with all individual files).

Be sure that you have made at least one reliable **backup** of all your old information in the old Radar program. Please read the topic **Automatic backup** in the Freenotes. See also the instruction video **Making backups** or read the chapters on making a backup in the Repertory utilities manual.

Import Radar Clipboards (Catalogue of cases)

From the menu **Tools**, select **Import**, then select the following option:

- **List of analyses from another program**

If in Radar you have used the default method of saving the clipboards, then in the next window at the option ‘*Choose the program that contains your analysis*’ you need to select:

- **Radar 10.5 catalogue of cases***

(* Not directly possible on Mac: see instruction)

To import you also need to select the location where the catalogue of cases is stored and select the “Catalogue name”.

*(You will only be asked to select 1 of the 4 files of a catalogue: the **xxxxx.da1** file. By default, it is the file: **repert.da1**).*

Then click on **Start import** to continue.




After the import of the catalogue of old Radar cases is ready, you can click on the button **Show my catalogues** to see the imported files.

Note: If you were using **several catalogues** of cases in the old Radar, you will need to repeat the import for each single catalogue.

When you import more than one catalogue of old cases, then a separate group (new catalogue) for each catalogue will be created in RadarOpus.

*Read more under **Recall and Save analysis in RadarOpus** for further explanation.*

Note: For most users just keep the option “**Import rubrics in national language if present**” selected.



Import Radar clipboards (separate Cas files)

Go from the menu **Tools** to **Import**, then select

- **List of analyses from another program**

In the next window select to import:

- **Radar single cases**

You also need to select the location of the **folder** where the single case files are located.

Advised settings:

Keep the import location set to: **Default folder**.

See the settings as shown on the screen:



If you want the cases stored in a separate folder in RadarOpus, then select: **Create a new folder**.

Note: Stored as external files: This option is **not** advised for the average user.

If the imported single cases should not really be imported into the new RadarOpus database, they again can be **stored as external, single files** (for each analysis a new separate file will be created on the hard disk). With this method, the separate files are converted to RadarOpus but remain single files on the hard disk.

After the import is completed you can click on the button **Show my folder** to Recall one of the imported cases in RadarOpus clipboards.

Read more below in the section **Recall and save analyses in RadarOpus**.



Import only One single Radar Case


If you do not want to import all your old Radar cases (clipboards) but only want to import one single case from time to time, this is also possible.

First export all your Radar cases to single CAS files. In your old Radar program go to **File**, then **Export**, then **Export all cases as single files (Cas files)**.

Inside the folder **C:\Radar\Cases** on your hard disk you will find a separate folder named "Reper" which will contain all your cases.


Note: If you use several catalogues of cases, you will need to export each one of them separately.

If you installed RadarOpus on a new computer copy this folder of CAS files to this new computer.



How to import a single case in RadarOpus:

- **Right** mouse click on a clipboard,
- then select **Recall analysis**.
- Here you can select to import a CAS file



Convert all your CAS files to RadarOpus XML files

If you prefer to keep on using individual files on your hard drive for each repertorisation (Analysis), you can convert all your CAS files to XML files.

- XML files is the new file type of RadarOpus which can be used to exchange (send) analyses for example by email to a colleague.
- XML files are also used by the Synthesis App. You can open or import such Synthesis App analysis also in your RadarOpus.

You could create for example on your local server a folder where you store these XML files. This will allow you to open such analyses from any workstation on which you have RadarOpus running

Go from the menu **Tools** to **Import**, then select

- **List of analyses from another program**


In the next window select to import:

- **Radar single cases**

You also need to select the location of the **folder** where the single case files are located.

Also select:

- **Save imported analyses in individual files** (instead of saving them in my RadarOpus database)




Attention: all your individual CAS files are converted to XML files, and you can find these in the folder: **\RadarOpus\Data\User\Singlecases**

You can move these converted XML file to any other location you prefer.

In RadarOpus you can recall and save such external XML files with the option:

- **Import an analysis**
- And
- **Export an analysis**



Recall (and Save) analyses in RadarOpus

After you have imported the old Radar cases you can go to the **Recall Analysis** window via the button **Show my folders**.

See video: [Save the analyses separately](#)


Note: “Cases” is now called “Analysis”

In RadarOpus a set of saved clipboards is called an ‘**Analysis**’.

The new RadarOpus clipboards can contain repertory symptoms and also symptoms taken from a Materia Medica or proving books. You can also take any other part of the RadarOpus Library into a clipboard. Anything taken into a clipboard can act as a link which will bring you back to the original location you took it from.

This **Recall Analysis** window (and also the **Save Analysis** window) can also be opened via:

- a **Right mouse click** on one of the clipboards.
- or via the **Analysis button** in the main toolbar.




Note for Mac users: see also in the chapter *Tips for Macintosh users* how to set the Function keys and the Right mouse click on a Mac.

All your old imported **Radar cases** can be found in this **Recall Analysis** window. There are several options in this window:

If you imported a “catalogue of cases” or “single Cas files” they can be found in this **Recall a case** window, via the option:

- **Recall an analysis from a folder**




Default folder

If you imported (and created) **several folders of cases** you can also find them here in this window via the **Folder selection button**. The default location for the RadarOpus Analysis is called **Default folder** of analyses. The by default imported set of **old Radar cases** (catalogue system) has the name **Reper**.

Save an Analysis

See the chapter *Save an Analysis* for more information on how to save your clipboards in RadarOpus.

In RadarOpus the clipboards are called ‘**Analysis**’



In RadarOpus two methods to save the Analysis

There are two methods of saving an analysis (clipboards) in RadarOpus. Please also see the RadarOpus manual for more information.

1. Much like in the old Radar this is called: **Save Analysis to a folder.**

See video: [Save the analyses separately](#)

2. Or saved in the: **Patient file.**

See video: [Save analyses to the patient file](#)

For more information see in the manual the chapters about *The Symptom Clipboards* and in *The Patient file*.


See video: [The Workflow of the patient file](#)

Import Winchip patient information from Radar


If you have used the **Winchip** patient file in Radar, you can import your patients' catalogue into RadarOpus*.

Be sure you have made at least one reliable backup of the folder **c:\radar\winchip\filepati** (with all sub folders) on a secondary media like an external hard-drive or memory stick.

Go to the menu **Tools**, then to **Convert/Import from other programs**, then **WinChip 3.x**



Select via the button **Browse** the location of the folder** **C:\Radar\Winchip**



Also select what you want to import. The average user would probably only like to import the actual Patient-database located in **C:\Radar\Winchip\Filepati**

(* **Note:** Switch from Windows to Mac:

If you switch from an old Radar **Windows** version to a **Mac** version of RadarOpus, your catalogue must be imported in a RadarOpus Windows version first. Then a RadarOpus backup should be made. This RadarOpus backup you can then import into RadarOpus. Your RadarOpus representative will assist you.

(** **Note:** You do not really need to select a folder with the exact name \Winchip, as long as this folder contains the sub-folder \Filepati (and optionally also \Uti and \Pck).

If you have only made a backup of the actual patient-database folder, then you need to place it in a folder with the name \Filepati

The imported patients will be found in the **Table Of Contents of Patients** (the list of Patients)

Import personal repertory additions from Radar

The option to import this information will be available at a later stage. Be sure you made and keep a reliable backup of folder **c:\radar\log**

Import personal keynotes, freenotes, EH-documents

The option to import this information will be available at a later stage. Be sure you have made a reliable backup of these files.

Important: Make regularly a RadarOpus backup

For more information please see the chapter Making Backups in the RadarOpus manual.